

Leo "S"


PANTOGRAPH CONSTRUCTION

Split side rails



Key features:

- Split side rails
- Easily removable bed ends
- Modern, pantograph construction
- Wide range of bed ends
- Wide height adjustment (lowest position 30cm)
- Bed ends independent from construction
- O Max. working load 215 kg
- 2 years guarantee







HEIGHT ADJUSTMENT

While designing Leo we paid careful attention to bed height adjustment - it is almost 50cm. Lowest position -30cm provides comfortable usage of bed and also reduces risk of falling. Maximum height - almost 80cm-remarkably ease care of lying patients for medical staff.







BED ENDS

are made of laminated board and solid beech wood.

Due to bed construction we can offer detachable foot/head ends.

The newest solution – Split Side Rail System. New wooden split side rails are divided 40/60 and protect the patient up to 400 mm of height. According to patients needs there is possibility of lifting part of the side rail from the head side or from the foot side. Split Side rails type is a solution which protects the patient during the operation and provides comfort and feeling of patient's safety for medical staff. New Split side rails are manufactured according to current standard EN-60601-2-52.







In our offer we have both plastic and metal castors, with various diameters.

For 50mm we provide possibility to order hidden castors. Mobility of the bed is increased by additional options, like individual (standard), axial or

Basic dimensions	
Overall width x lenght	
Mattress platform width x lenght	



adjustment



section

adjustment



Auto-Contour

function



angle 15º

central brake.



Trendeleburg position's angle 17

Height adjustment 50cm

Adjust design to individual needs

Wide range of bed ends, shapes and great choice of finishings



Bicoloured wooden bed end



L-04 Carved wooden bed end

127 x 207 110 x 200



L-05 Solid wooden bed end





Standard colours





Non-standard colours





