

LEOL Rehabilitation bed on pantograph construction Full lenght side rails



LEO L is a modern rehabilitation bed offering a wide range of conveniences for both the user and the caregivers. It meets the highest world's standards and provides full support for disable patients and people who need long term care. Leo "L" is basic model in Leo family and what distinguishes it from other variants is full lenght side rails.

- Full lenght side rails
- Modern, pantograph construction
- Wide range of height adjustment: lowest mattress platform position ~31cm
- Bed ends independent from construction
- Wide range of bed ends, shapes and great choice of finishings
- Intuitive and easy to use 10 buttons handset
- Max. working load: 215kg
- 2 years guarantee



REHA-BED Sp. z o.o. ul. Spacerowa 1, 41-253 Czeladź





) biuro@rehabed.com.pl





Leo L Modern construction and attractive design

WOODEN FULL LENGHT SIDE RAILS are the standard equipment of all nursing care beds. They protect the patient up to 385 mm over the mattress platform. For additional protection there is available a set of Side Rails Height Extension, which increases the height of side rails up to 530 mm. This solution is very helpful when the use of pressure mattresses is needed. Full-length wooden side rails are manufactured according to current standard EN-60601-2-52

HEIGHT ADJUSTMENT

While designing the Leo L, we paid carefull attention to the range of height adjustment of the bed frame. For the Leo L rehabilitation bed, it is nearly 50cm. The lowest position of mattress platform is ~31cm (varies depending on the

diameter of the wheels) It allows the patient to get out of bed comfortably and significantly reduces the risk of falling. The highest position, nearly 80cm. remarkably ease care of lying patients for medical staff.



HOLDERS FOR ADDITIONAL EQUIPMENT

Leo beds are standardly equipped with special sockets, located both on left and right side of the head section of bed frame. Those holders enable assembly of additional equipment - lifting and transfusion pole.



HANDSET provides comfortable one-hand usage and fast access to all bed functions. It is also possible to block individual functions.

CASTORS

We offer plastic and metal wheels with a choice of diameter. For wheels with a diameter of 50mm, there is an option to order a

concealed variant. Additional options that enhance the mobility of the bed include individual wheel locking (standard), axial locking, and central locking.

Auto-Contour

function



215ka

Max. working load



Basic dimensions		
Overall width x lenght	105 x 206cm	
Mattress platform w. x l.	90 x 200cm	
		Backrest section

Adjust to individual needs

BED ENDS

Bed ends made of laminated board and solid beech wood.

Standard wood colours





Tigh and calf

section adjustment

adiustment

- Trendeleburg 15°
- Reverse Height Trendeleburg 179

adjustment 31.5-78.5cm

- L-06 Rounded full bed end L-05 Full bed end
- L-04 Curved bed end
- L-03 Two-tone wooden bed end
- L-02 Wooden bed end with wooden slat

L-01 Wooden bed end with chrome tube

WIZUALIZACJA ŁÓŻKA REHABILITACYJNEGO JEST GRAFIKĄ POGLĄDOWĄ, SZCZEGÓŁY WYKOŃCZENIA ORAZ WYPOSAŻENIA ZAWSZE SĄ INDYWIDUALNIE USTALANE Z DZIŁAEM HANDLOWYM. RZECZYWISTE KOLORY MOGĄ RÓŻNIĆ SIĘ OD PRZEDSTAWIONYCH NA ZDJĘCIACH.

REHA-BED Sp. z o.o. ul. Spacerowa 1, 41-253 Czeladź



biuro@rehabed.com.pl

